



CERTIFICATION SHEET

Umbilical Pressure Leak Test

| Umbilical Serial # | System Component Part Description /Identification |
|--------------------|---|
| | |
| | |
| | |

| Purpose of Tests (check test box) |
|---|
| <input type="checkbox"/> Conventional Umbilical Routine Pressure Test , as well as visual inspection of umbilical for cuts, slices, bulges, fitting slippage, excessive wear, signs of kinking, and any other signs of damage or deformity, 600 psig for 10 min. |
| <input type="checkbox"/> XLDS Umbilicals Routine Pressure Test, as well as visual inspection of umbilical for cuts, slices, bulges, fitting slippage, excessive wear, signs of kinking, and any other signs of damage or deformity, 800 psig for 10 min. |
| |
| |

| Maximum System Design Pressure | Maximum Umbilical Rated Working Pressure |
|--------------------------------|--|
| psig | psig |

| Test Pressure Used | Test time Used | Test Media Used |
|--------------------|----------------|--|
| psig | | <input type="checkbox"/> Water <input type="checkbox"/> Divers Air <input type="checkbox"/> N2 <input type="checkbox"/> Other |

| Test 1 |
|---|
| Start Time _____ Start Pressure _____ Stop Time _____ Stop Pressure _____ <input type="checkbox"/> Pass <input type="checkbox"/> Fail |

| |
|------------------|
| Remarks (if any) |
| |
| |

| Print Name | Signature | Date |
|------------|-----------|------|
| | | |